




# CHROMIUM, HEXAVALENT

(0 to 0.60 mg/L Cr<sup>6+</sup>) Method 8023

For water and wastewater

1,5-Diphenylcarbohydrazide Method\* USEPA accepted for wastewater analyses\*\*



1. Enter the stored program number for hexavalent chromium (Cr<sup>6+</sup>)- powder pillows.

Press: **PRGM**

The display will show:

**PRGM ?**



## CHROMIUM, HEXAVALENT



2. Press: **13 ENTER**

The display will show  
**mg/L, Cr6** and the  
**ZERO** icon.

*Note: For alternate forms  
( $\text{CrO}_4$ ,  $\text{Cr}_2\text{O}_7$ ), press the  
**CONC** key.*



## CHROMIUM, HEXAVALENT



**3.** Fill a sample cell with 10 mL of sample.



## CHROMIUM, HEXAVALENT




**4.** Add the contents of one ChromaVer 3 Reagent Powder Pillow to the cell (the prepared sample). Cap the cell and invert several times to mix.

*Note: A purple color will form if  $\text{Cr}^{6+}$  is present.*



## CHROMIUM, HEXAVALENT



5. Press:

**TIMER ENTER**

A five-minute reaction period will begin.



## CHROMIUM, HEXAVALENT




**6.** Fill another sample cell with 10 mL of sample (the blank).

*Note: For turbid samples, add the contents of one Acid Reagent Powder Pillow. This ensures turbidity dissolved by the acid in the ChromaVer 3 Chromium Reagent is also dissolved in the blank.*




## CHROMIUM, HEXAVALENT



7. When the timer beeps, place the blank into the cell holder. Tightly cover the sample cell with the instrument cap.



## CHROMIUM, HEXAVALENT




8. Press: **ZERO**

The cursor will move to the right, then the display will show:

**0.00 mg/L Cr6**



## CHROMIUM, HEXAVALENT



9. Place the prepared sample into the cell holder. Tightly cover the sample cell with the instrument cap.

# CHROMIUM, HEXAVALENT



## 10. Press: **READ**

The cursor will move to the right, then the result in mg/L hexavalent chromium will be displayed.

*Note: Standard Adjust may be performed using a prepared standard (see Standard Adjust in Section 1).*

