




CHEMICAL OXYGEN DEMAND, HR

Colorimetric Determination, 0 to 1,500 and 0 to 15,000 mg/L COD



1. Enter the stored program number for chemical oxygen demand (COD), low range.

Press: **PRGM**

The display will show:

PRGM ?



CHEMICAL OXYGEN DEMAND, HR



2. Press: **17 ENTER**

The display will show
mg/L, COD and the
ZERO icon.

Note: For alternate form (O₂), press the
CONC key.



CHEMICAL OXYGEN DEMAND, HR



3. Insert the COD/TNT Adapter into the cell holder by rotating the adapter until it drops into place. Then push down to fully insert it.

Note: For increased performance, a diffuser band covers the light path holes on the adapter. Do not remove the diffuser band.



CHEMICAL OXYGEN DEMAND, HR



4. Clean the outside of the blank with a towel.

Note: Wiping with a damp towel, followed by a dry one, will remove fingerprints or other marks.



CHEMICAL OXYGEN DEMAND, HR




5. Place the blank in the adapter.

Push straight down on the top of the vial until it seats solidly into the adapter.

Note: Do not move the vial from side to side as this can cause errors.



CHEMICAL OXYGEN DEMAND, HR




6. Tightly cover the vial with the instrument cap.

Note: *The blank is stable when stored in the dark. See Blanks for Colorimetric Determination following these procedures.*



CHEMICAL OXYGEN DEMAND, HR



7. Press: **ZERO**

The cursor will move to the right, then the display will show:

0 mg/L COD



CHEMICAL OXYGEN DEMAND, HR



8. Clean the outside of the sample vial with a towel.



CHEMICAL OXYGEN DEMAND, HR




9. Place the sample vial in the adapter.

Push straight down on the top of the vial until it seats solidly into the adapter.

Note: Do not move the vial from side to side as this can cause errors.




CHEMICAL OXYGEN DEMAND, HR



10. Tightly cover the vial with the instrument cap.



CHEMICAL OXYGEN DEMAND, HR



11. Press: **READ**

The cursor will move to the right, then the result in mg/L COD will be displayed.

Note: When using High Range Plus COD Digestion Reagent Vials, multiply the reading by 10.

Note: For most accurate results with samples near 1,500 or 15,000 mg/L COD, repeat the analysis with a diluted sample.